

# CTS-10988

2690/2695 Performance Maintenance Kit  
Comparable to: Waters® WAT270944



## Technical Data Sheet

### 1.0 Description

This PM kit contains all the parts necessary to keep Waters® 2690, 2690D, 2695, 2695D and Alliance® HPLC systems working at peak performance. It is comparable to Waters® part number WAT270944.

### 2.0 Specifications

#### 2.1 Components

Qty	Part Description
2	Sapphire Plunger (CTS-10476)
1	Seal Wash Plunger Seal Kit, 2/pk (CTS-10560)
1	Head Plunger Seal Kit, 2/pk (CTS-10566)
1	Wash Tube Seal Kit, 4/pk (CTS-10597)
1	Sparge Diffuser (CTS-10984)
1	Face Seal, 4/pk (CTS-10989)
4	Solvent Reservoir Filter, 10µm (CTS-5531)
1	250µL Syringe (CTS-3109)
1	2µm Precolumn Filter Insert (CTS-10919)
2	Seal Wash Face Seal (CTS-10532)
1	PP Check Valve Cartridge, 2/pk (CTS-10600)
1	Seal Pack Rebuild Kit (CTS-10839)

#### 2.2 Wetted Materials

18-8 SS  
302 SS  
316 SS  
316 SST



Borosilicate Glass  
Delrin®  
Elgiloy Spring  
HDPE  
MP35N SS  
PCTFE  
PEEK™  
PTFE  
Ruby  
Sapphire

### 3.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.

### 4.0 Warranty

We view your business with us as a partnership based on mutual trust.

We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at [sales@sciencix.com](mailto:sales@sciencix.com) to report the issue. We will assist you with a resolution in a timely fashion.