

CTS-5531-05

Solvent Bottle Filter, 5/pk
Comparable to: Waters® 289003536



Technical Data Sheet

1.0 Description

These solvent bottle filters are designed to be compatible with 510, 515, 600, 610, 1515, 1525, 1525 Micro, 1525EF, 2535, 2545, 2555, 2690, 2690D, 2695, 2695D, 2700, ACQUITY UPLC® BSM systems. They are comparable to Waters® part number 289003536.

2.0 Wetted Materials Specifications

316 SS

3.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.

Please note: A short piece of 1/8" ID tubing, such as CTS-6809 (Mobile Phase Transfer Line, TFE, 3/16" OD x 1/8" ID x 25'), is needed to attach the CTS-5531 solvent filter to standard 1/8" OD solvent tubing.



4.0 Warranty

We view your business with us as a partnership based on mutual trust.

We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at sales@sciencix.com to report the issue. We will assist you with a resolution in a timely fashion.

Please note: A short piece of 1/8" ID tubing, such as 11-1112 (Tubing FEP 3/16" OD x 1/8" ID x 10ft), is needed to attach the CTS-5531 solvent filter to standard 1/8" OD solvent tubing.