

CTS-10919-5

2µm Precolumn Filter Insert, 5/pk
Comparable to: Waters®
WAT05139



Technical Data Sheet

1.0 Description

These precolumn filters are designed to be compatible with Waters® 600, 717, 717 Plus, 2690, 2690D, 2695, 2695D, 2790, 2795, 2796, Alliance®, CapLC and Cap LC XE HPLC systems. They are comparable to Waters® part number WAT05139.

2.0 Materials Specifications

316 SS

3.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.



4.0 Warranty

We view your business with us as a partnership based on mutual trust.

We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at sales@sciencix.com to report the issue. We will assist you with a resolution in a timely fashion.