

CTS-10527

Precolumn Filter Assembly
Comparable to: Waters®
WAT035190
Old# WAT084561



Technical Data Sheet

1.0 Description

This precolumn filter assembly is designed to be compatible with Waters® 600, 717, 2690, 2690D, 2695, 2695D, 2790, 2795, 2796 and Alliance® HPLC systems. It is comparable to Waters® part number WAT035190. It replaces old part number WAT084561.

2.0 Materials Specifications

316 SS
316L SS

3.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.



4.0 Warranty

We view your business with us as a partnership based on mutual trust.

We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at sales@sciencix.com to report the issue. We will assist you with a resolution in a timely fashion.