

# CTS-11434

Sapphire Plunger  
Comparable to: Shimadzu®  
228-35010-91



## Technical Data Sheet

---

### 1.0 Description

This plunger is designed to be compatible with Shimadzu i-Series, LC-2010 A/C HT, LC-HT, SIL-10ADvp, SIL-20A/AC, SIL-20ACXR, SIL-20ACHT, SIL-30AC, SIL-30ACMP, SIL-HTa and SIL-HTc HPLC systems. It is comparable to Shimadzu part number 228-35010-91.

### 2.0 Materials Specifications

Sapphire  
316 SS

### 3.0 Replacement Frequency

Replacement may be indicated by a leak at the pump head. The sapphire plunger should be inspected with a magnifier whenever the pump head is disassembled to replace the plunger seals. If any scratches, a worn area or flat area are seen then the plunger should be replaced.

### 4.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.



### 5.0 Warranty

We view your business with us as a partnership based on mutual trust.

We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at [sales@sciencix.com](mailto:sales@sciencix.com) to report the issue. We will assist you with a resolution in a timely fashion.