

# CTS-21936

Foam Swab, Anti-Static, 50/pk  
Comparable to: Sciex® 1017396



## Technical Data Sheet

---

### 1.0 Description

These swabs are designed to be compatible with 2000, 3000, 3200, 3200QTRAP 4000, 5000, 4000QTRAP, QTRAP4500, QTRAP, 4600, 5500+, 5600, 6600, QTRAP6500 systems. They are comparable to Sciencix part number 1017396.

### 2.0 Wetted Materials Specifications

Anti-Static Foam

### 3.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper usage procedure.



### 4.0 Warranty

We view your business with us as a partnership based on mutual trust.

We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at [sales@sciencix.com](mailto:sales@sciencix.com) to report the issue. We will assist you with a resolution in a timely fashion.

**Please note:** A short piece of 1/8" ID tubing, such as 11-1112 (Tubing FEP 3/16" OD x 1/8" ID x 10ft), is needed to attach the CTS-5531 solvent filter to standard 1/8" OD solvent tubing.