

CTS-21456

Performance Maintenance (PM) Kit VF-P10, HPG-R
Comparable to: Thermo™ Dionex™ 6040.1956B



Technical Data Sheet

1.0 Description

This kit contains all the parts necessary to keep Thermo™ Dionex™ VF-P10 and HPG-R systems working at peak performance. It is comparable to part number 6040.1956B.

2.0 Specifications

2.1 Components

| Qty | Part Description |
|---------|---|
| 2 | Valve Cartridge, Zirconia (CTS-21342) |
| 0.50ft | Pharmed Tubing 1/16" ID x 1/8" OD (51108) |
| 5 | 1/16" Double Barb Fitting - Clear PP (CTS-21443) |
| 1 | 1/16" Double Barb Fitting - 90 Clear PP (CTS-21442) |
| 3 | J-Style Clip Mount, Clear (CTS-21455) |
| 4 | Bearing Ring (CTS-21411) |
| 4 | Piston Seal Assembly (CTS-21318) |
| 1 | Ring Filter (CTS-21299) |
| 2 | 9mm X 1.5mm PTFE O-Ring (CTS-21412) |
| 2 | 32mm X 1.5mm PTFE O-Ring (CTS-21413) |
| 1 | Purge Cap (CTS-21294) |
| 4.92 ft | 1.3mm ID X 2.8mm OD Silicone Tubing (51109) |
| 4 | Disk Filter (CTS-21298) |
| 4 | Filter Holder Assembly (CTS-21457) |
| 5 | Puritan 892-PC DBL PurSwab (51123) |



2.2 Materials Specifications

PTFE
Zirconia
Clear PP
Clear PVC
PEEK
Titanium
Silicone
Ceramic
UHMW-PE
316 SS

3.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.

4.0 Warranty

We view your business with us as a partnership based on mutual trust. We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at sales@sciencix.com to report the issue. We will assist you with a resolution in a timely fashion.