

# CTS-21501

Maintenance Kit, w/o Viper™ Capillary,  
Compatible with Thermo™/Dionex™ WPS 3000 SL/RS,  
Comparable to Thermo™/Dionex™ # 6840.1951C, Old #  
6840.1951B



## Technical Data Sheet

### 1.0 Description

This kit contains all the parts necessary to keep Thermo™ Dionex™ WPS 3000 SL/RS systems working at peak performance. It is comparable to Thermo™/Dionex™ OEM # 6840.1951C, Old # 6840.1951B.

### 2.0 Specifications

#### 2.1 Components

| Qty | Part Description                               |
|-----|------------------------------------------------|
| 1   | Service Kit, Needle Seat (CTS-21502)           |
| 1   | Seal Washer, Syringe, WPS 5/pk (CTS-21571)     |
| 1   | Sample Needle, Viper™ (CTS-21505)              |
| 1   | Peristaltic Pump Cartridge, WPS/AC (CTS-21406) |
| 1   | 100µL Syringe for WPS-3000 Series (11-3332)    |
| 1   | Rotor Seal Kit (CTS-21507)                     |



### 3.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.

### 4.0 Warranty

We view your business with us as a partnership based on mutual trust. We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at [sales@sciencix.com](mailto:sales@sciencix.com) to report the issue. We will assist you with a resolution in a timely fashion.

### 2.2 Materials Specifications

PEEK  
316 SS  
PTFE  
Titanium Nitride  
White Derlin®  
PharMed® BPT PP-TPE  
Borosilicate Glass  
303 SS  
304 SS  
PCTFE  
VespeI®