

CTS-22054

Performance Maintenance Kit

Compatible with Thermo™/Dionex™ VH-P1 and HPG-SD,
Comparable to Thermo™/Dionex™ # 6044.1956A



Technical Data Sheet

1.0 Description

This kit contains all the parts necessary to keep Thermo™ Dionex™ VH-P1 and HPG-SD systems working at peak performance. It is comparable to part number 6044.1956A.

2.0 Specifications

2.1 Components

Qty	Part Description
6	Filter, 10 Micron, Titanium (CTS-22081)
5	Filter Holder Assy (CTS-21457)
2	9mm X 1.5mm-PTFE O-Ring (CTS-21412)
2	O-Ring, PTFE, 45mm ID x 1.5mm CS (CTS-22390)
2	O-Ring, PTFE, 65mm ID x 1.5mm CS (CTS-22391]
4	Piston Seal, Analyt, RP without Flange, Biocomp (CTS-22038)
4	Bearing Ring (CTS-21411)
1	Pharmed Tubing, 1/16" ID x 1/8" OD x 32.5" Long (CTS-22393)
4	Tubing Fitting, 1/16" Barb to M5 x 0.8 mm Thread, PP (CTS-22372)
1	1/16" Double Barb Fitting-90 Clear PP (CTS-21442)
1	1/16 Double Barb Fitting - Clear PP (CTS-21443)
2	J-Style Clip Mount, Clear, Adh Mnt (CTS-21455)
1	Silicone Tubing, 1.3mm ID x 2.8mm OD x 60" Long (CTS-22394)
8	Screw, M5 x 0.8 x 16mm, SST, Pan Head, Torx (CTS-22392)
25	Lint-Free Swab (51123)



2.2 Materials Specifications

Titanium
PEEK
PTFE
UHMW-PE
PharMed®
Polypropylene
PVC
Silicone
18-8 SS
Cotton
Fine Paper

3.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.

4.0 Warranty

We view your business with us as a partnership based on mutual trust. We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at sales@sciencix.com to report the issue. We will assist you with a resolution in a timely fashion.