

# CTS-21692

Acquity® ISM Performance Maintenance Kit  
Comparable to: Waters® 201000286



## Technical Data Sheet

### 1.0 Description

This kit contains all the parts necessary to keep Waters Acquity® ISM systems working at peak performance. It is comparable to Waters® part number 201000286.

### 2.0 Specifications

#### 2.1 Components

Qty	Part Description
1	SS Filter Frit Replacement Cartridge, 2/pk (CTS-11415)
1	Acquity H-Class Filter Assy, 22uL (CTS-21196)
1	UPLC Check Valve Assy, Spring Loaded (CTS-21100)
1	Plunger, 2/pk (CTS-11093-02)
1	O-Ring, PTFE (CTS-11092)
1	Seal Wash, 2/pk (11-4181)
1	Head Plunger Seal, 2/pk (CTS-11119-02)
1	Tubing, Teflon (CTS-21138)
2	Solvent Filter, 10 µm (CTS-21256)
1	Filter, Air, Door (CTS-11370)
1	Check Valve, Double Ball & Seat (CTS-A11095)

#### 2.2 Materials Specifications

316 SS

Ruby

Zarconia



Teflon  
PEEK™  
Polyurethane Foam  
PTFE  
Synthetic Sapphire  
UHMW-PE  
MP35N  
Alumina

### 3.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.

### 4.0 Warranty

We view your business with us as a partnership based on mutual trust.

We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at [sales@sciencix.com](mailto:sales@sciencix.com) to report the issue. We will assist you with a resolution in a timely fashion.