

# CTS-11141

Acquity® SM Performance Maintenance Kit

\*Injector Valve Cartridge (700005236) not included

\*Rotor Seal & Pod Repair Kit included

Comparable to: Waters® 201000174



## Technical Data Sheet

### 1.0 Description

This kit contains all the parts necessary to keep Waters Acquity and Acquity UPLC® Sample Manager systems working at peak performance. It is comparable to Waters part number 201000174.

### 2.0 Specifications

#### 2.1 Components

Qty	Part Description
2	2.5mL Wash Syringe (11-3595)
1	Rotor Seal & Pod Repair Kit, Vespel® 3-Groove for Acquity Everflow Injector Pod/Cartridge (700002765) (11-4134)
1	Needle Seal O-Ring, 002 Kalrez® (CTS-11149)
1	Filter, Air, Chassis (CTS-11150)
1	Tube Assy, Inject Outlet (CTS-A11142)
1	Tube Assy, System Outlet, SDS (CTS-A11145)
1	Sample Loop, 10µL (CTS-A11147)
1	Needle Assembly, PEEK™ (CTS-A11151)
1	100µL Sample Metering Syringe, HP (CTS-A11167)

#### 2.2 Materials Specifications

303 SS

304 SS

316 SS

Borosilicate Glass

Glass

Kalrez



Neoprene

PEEK

Polyurethane Foam

PTFE

UHMWPE

Vespel

### 3.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.

### 4.0 Warranty

We view your business with us as a partnership based on mutual trust.

We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at [sales@sciencix.com](mailto:sales@sciencix.com) to report the issue. We will assist you with a resolution in a timely fashion.