

CTS-11380

2690/2695 Performance Maintenance Kit w/Seal Pack



Technical Data Sheet

1.0 Description

This PM kit contains all the parts necessary to keep Waters 2690, 2690D, 2695, 2695D and Alliance HPLC systems working at peak performance.

2.0 Specifications

2.1 Components

Qty	Part Description
2	Sapphire Plunger (CTS-10476)
1	Seal Wash Plunger Seal Kit, 2/pk (CTS-10560)
1	Head Plunger Seal Kit, 2/pk (CTS-10566)
1	Wash Tube Seal Kit, 4/pk (CTS-10597)
1	Seal Pack Assembly Kit (CTS-11222)
1	2µm Precolumn Filter Insert (CTS-10919)
1	Sparge Diffuser (CTS-10984)
1	Face Seal, 4/pk (CTS-10989)
1	250µL Syringe (CTS-3109)
4	Solvent Reservoir Filter, 10µm (CTS-5531)
1	PP Check Valve Cartridge, 2/pk (CTS-10600)



Borosilicate Glass

Delrin®

Elgiloy Spring

HDPE

MP35N SS

PCTFE

PEEK™

PTFE

Ruby

Sapphire

3.0 Installation Information

Sciencix parts are designed to be compatible and interchangeable with OEM parts. Please see instrument service manual for proper removal and installation procedure.

4.0 Warranty

We view your business with us as a partnership based on mutual trust.

We extend a Lifetime Warranty against any manufacturing defects on all Sciencix HPLC repair parts. If a defect is identified, please contact us at 800.682.6480 or by email at sales@sciencix.com to report the issue. We will assist you with a resolution in a timely fashion.

2.2 Materials

18-8 SS

302 SS

316 SS

316 SST